

Bacterial Strain Preservation

BacterialStrainPreservation.doc
Last Updated: 080723 by Janet Bautista

INTRODUCTION

The method described is to preserve bacteria strains of *Pseudomonas aureofaciens* and *Pseudomonas chlororaphis* in 15% glycerol solution for long-time storage in a -80°C freezer. The purpose of storing strains of bacteria in a -80°C freezer is to have archives of newly acquired or created strains.

SUPPLIES NEEDED

- Fresh liquid bacteria culture
 - o Grown overnight (at least eight hours)
- Sterile 60% glycerol solution (autoclave)
- micropipette
- 2 sterile pipette tips
- sterile micro centrifuge
- liquid nitrogen
- -80°C freezer

ROUTINE PROCEDURE

Making 60% Glycerol Solution

- 1) The solution is stored in the refrigerator.
- 2) If solution is low, make more by using a purple screw cap 100 ml media stock bottle to store the solution.
- 3) Using a graduated cylinder, measure 30 ml of Glycerol. Pour into the media bottle.
- 4) Rinse the graduated cylinder with DI water and then fill to 20 ml of DI water. Pour into the media bottle. The total volume of solution should be 50 ml.
- 5) Vortex or just swish the solution well and have it autoclaved.
- 6) Store solution in the refrigerator.

Making Fresh Liquid Bacteria Culture (Inoculation)

P. aureofaciens + *P. Chlororaphis*

- 1) Use flamed loop to transfer a single colony of *P. aureofaciens* from plate #2 or #3 to 5 ml of nutrient broth in a sterilized culture tube.
- 2) Leave on shaker overnight

Autoclaving

- 1) Autoclave as part of a larger batch for efficiency at 121°C for 30 minutes.
 - a. Located in Johnson Hall room 227
 - b. Autoclave should be set to 'liquid 60'.
- 2) Supplies to autoclave include:
 - a. Micro centrifuge tubes
 - b. Pipette tips

- c. 60% glycerol solution
- d. Include other supplies related to method for efficiency (ie: Petri plates, etc)

Making the Bacteria Archive

- 1) Label a sterile micro centrifuge tube with the date and name of bacteria
- 2) Pipette using a sterile pipette tip 0.1 ml of **sterile 60% glycerol solution** into the tube.
- 3) Using a new pipette tip, add 0.3 ml of the bacterial culture (frozen stock will be 15% glycerol). Use a new pipette tip for each bacteria strain.
- 4) Cap the tube and shake to mix the bacteria culture and the glycerol solution together.
- 5) Using a cooler filled with liquid nitrogen, immediately place the tubes inside to flash freeze the bacteria stock.
- 6) Bring cooler to Oceanography building (Marine Sciences) room 207 where the -80°C freezer is located and store tube in the freezer.
- 7) Remember to bring back the cooler.

Retrieving Bacteria Culture from Freezer Stock

- 1) In the Oceanography building (Marine Sciences) room 207/257, the stock tubes of the bacteria strains are kept in a -80°C freezer.
- 2) Label a Petri plate with bacteria species and date.
- 3) Use a sterilized toothpick to streak cells (plate #1).
- 4) Incubate the plate in the bacterial incubator located in room 303B.
- 5) Two to four days later, transfer a single colony by flamed loop from plate #1 to new plate and streak (plate #2).
- 6) Two to four days later, transfer a single colony by flamed loop from plate #2 to new plate and streak (plate #3).
- 7) Use only plate #2 and plate #3 for inoculation.

TROUBLESHOOTING

- 1) If you mess up start over.

REFERENCES

- 1) *Long Term Storage of Bacterial Strains*, Nov. 14, 1990. Donis Keller Lab, C. Helms.
- 2) *Bacterial Glycerol Stocks*, University of North Carolina.
- 3) *Glycerol Stocks of Transformed Cells*, South African Structural Biology Initiative.